

## Demountable Beta Gas Cell 5cm or 10cm



The demountable Beta Gas Cell comes complete with a 10cm or 5cm Pyrex Cell body, cell mount, 2 Viton® O-rings or gaskets, 2 PTFE gaskets, and 2 septa, but without windows. (Windows and gaskets can be ordered separately). The window size is 25mm x 4mm. When assembled, the cell will be gas tight and operable at both reduced or positive pressures. Filling of the cell is conveniently handled with a gas tight syringe (ICL Part no. 0009-609). The exploded view shows the assembly and disassembly sequence of the cell.

0001-1469 PTFE Gaskets - 25mm dia, pk/4

0001-3729 Viton® Gaskets - 25mm dia, pk/4

0008-504 Demountable Beta Gas Cell- 10cm

0008-505 Demountable Beta Gas Cell- 5cm

0008-616 Septa Replacement Caps for Demountable Beta Gas Cell- pk/6

0008-4469 Split Retainer Rings for Demountable Beta Gas Cell- pk/2

0008-4452 Beta Gas Cell Holder Only



**JASCO France**  
17 Rue des Cerisiers  
91090 Lisses - France

Tél. 01 64 97 09 60 Fax: 01 85 09 90 76  
email: [jascofrance@jascofrance.fr](mailto:jascofrance@jascofrance.fr)