

Ideal for applications such as: Cement, Metals, Mining & Food Stuffs



## Product Features

- Polycarbonate safety guards
- Adjustable upper bolster
- Pressure release valve
- Adjustable pressure control valve
- Vacuum ports
- Gauges for low pressure applications (0-1, 0-2, & 0-5 ton options)



## Ordering Information:

GS15011 Atlas 15T Manual Hydraulic Press

GS25011 Atlas 25T Manual Hydraulic Press

### Spares and Consumables

GS15100 Replacement Seals and Gaskets kit for Manual Hydraulic Presses

GS15101 Hydraulic Oil (1litre)

### Low Tonnage Gauge Kits

GS15051 0-1 Ton Gauge Conversion Kit

GS15052 0-2 Ton Gauge Conversion Kit

GS15055 0-5 Ton Gauge Conversion Kit

## Specifications

Max height (at handle)	610mm
Width	310mm
Depth	190mm
Weight	50kg
Lower piston stroke	25.4mm
Upper bolster screw travel	89mm
Distance between pressing faces (max)	152mm
Distance between pressing faces (min)	38mm
Lower pressing face diameter	86mm
Upper pressing face diameter	32mm
Max width of sampling area (side to side)	143mm
Max depth of sampling area (front to back)	141mm